

OPTIMIZED APPLICATORS



Micro
Needle



Diamond
Peel



Aqua
Peel



LED Light
Therapy

AQUAPUREII FUNCTIONS INCLUDE

- Skin Rejuvenation

Deep Exfoliation & Extraction

Skin Texture

Transdermal Delivery
- Skin Whitening

Skin Boost

Acne & Anti-Bacteria

Skin Relaxation

SYSTEM SPECIFICATIONS

Input Voltage	100~240V, 50/60Hz
Panel	Touch panel
Dimension	365 x 523 x 460mm (WxDxH)
Weight	14kg

Entire Beauty Solution



EXPORT ONLY

EXPORT ONLY



AQUAPURE II

SMART GENUINE FACIAL CARE SYSTEM


Go Pure Be Youth

AQUAPURE II




DEEP EXFOLIATION BY AD PEEL TECHNOLOGY

AD Peel Technology produces channels on the skin to enhance delivery of solutions for maximizing the effects of microdermabrasion.



VACUUM ASSISTED MF(MICRO-FOCUSED) NEEDLE

New transdermal delivery technology with vacuum maximizes treatment comfort, efficiency and results.



New Upgraded AD Peel Technology

Advanced concept of AD(Aqua+Diamond) peel technology improves skin exfoliation and extraction in a shorter time, thus reducing treatment time.

Precise Solution Control Technology

Utilizes a stable adjustment valve to deliver equal distribution of all solutions being applied during all treatment circumstances.

Simultaneous Use with LED Light

Operator is able to treat two patients simultaneously.

User-friendly Smart Auto Reaction

The Smart Auto Reaction Interface detects and enables automatic recognition of all activated hand-pieces.

Higher Suction Power

Equipped with dual motor heads to produce better treatment results.

*2.5 times stronger suction power than previous version.

LIGHTPURE LED COLOR THERAPY


625nm

Cell activity
Skin firming anti-aging efficacy




520nm

Balance the water and oil



465nm

Inflammation
Bactericidal and anti-inflammatory



AQUASOLUTION



PEEL+

The first step of treatment is the process of non-irritating skin scaling which facilitates accelerated exfoliation, the removal of dead skin cells and other impurities.



REJUVE+

The final step of treatment promotes skin elasticity by inducing collagen and elastin regeneration and improving the inhibition of moisture loss by improving water retention.



SEBO+

The second step of treatment is for sebum dissolution, treatment of blackheads, and effective pore control. This assists the skin to balance between oily and moisture textures.



CLEANSYS

To provide a convenient cleaning agent specifically made for the needs of the system, CLEANSYS assists in easy maintenance procedures for the AQUAPURE.

CLINICAL CASES



Baseline

Immediately Post Procedure



Baseline

Immediately Post Procedure